





How to update the VSP for Twain drive

- Date : 2021.09.13
- Ver. : 1.0

1. Twain update


1-1. Close VSP Tray



1. Quit all kinds of VSP TRAY on task bar first.
✗ Close EzDent-i program.

1. Twain update

1-2. Run the VSP twain drive installer file



1. Run as administrator the VSP twain drive.
(VSP_Twain_driver_installer_1.0.x.xxx.exe)

※ Twain version



※ Driver installer storage(USB)
: Twain installation file is included in the USB


1. Twain update

1-3. Update VSP twain drive




1. Twain update

1-3. Update VSP twain drive



1. Twain update

1-3. Update VSP twain drive



1. Installation is complete. Click the Finish button.